

Rigid Bronchoscopy for Cyrobiopsy with Jet Ventilation

A combined approach using a flexible bronchoscope introduced through a rigid bronchoscope during cryosampling which provides both the flexibility to reach less accessible regions in the tracheobronchial tree and the ability to rapidly withdraw the flexible bronchoscope cryoprobe and attached frozen tissue, and re introduce the scope.

As the frozen tissue attached to the probe is large, cryoprobe and adhered tumor are removed together quickly; this causes the bronchoscope to move significantly, and often an assistant is needed to hold the bronchoscope permanently, to keep it in place for further introduction of the cryoprobe at the right spot! Use of Mediflex's BSS will assist here significantly.



"For several weeks I had the opportunity to test the Mediflex Bronchoscopic Stabilizer. The stabilizer is made of solid steel. At the same time, the place in contact with the bronchoscope is covered with soft Silicone Liners. The equipment is solid, you can trust it holding a bronchoscope. I used it for several classical bronchoscopy procedures clearing the trachea and placing a stent. The tool was particularly helpful during electromagnetic bronchoscopic navigation with cryo-biopsy of a lung tumor."

- Dr. Robert Dzedzic , Interventional Pulmonologist
Med. Univ. of Gdansk (MUG) Dept. of Thoracic Surgery