



Accuracy Begins With Stability

ELIMINATE AN ASSISTANT

Hands-Free static bronchoscope positioning eliminates the need for an assistant, and the associated fatigue and tremor

PHYSICIAN (PROCEDURE) FOCUS

Physicians can focus on the procedure, without having to manage the bronchoscope

REDUCE RADIATION EXPOSURE

Allows the physician and team to step away during imaging

GREATER ACCURACY (ENHANCED PRECISION)

Supports accuracy of diagnosis, intervention, and tissue collection with decreased sampling error

SHORTEN PROCEDURE

Maintains optimal bronchoscope positioning resulting in greater efficiency thus reduced procedure time